

Measuring Cylinders

Conical shaped cup with graduated markings on the side to allow easy and accurate measurement of volume of liquids

- Made of Borosilicate Glass
- Amber Stain Graduation
- With Spout and Solid Base



Item No.: 1612-01



Item No.: 1612-B



Item No.: 1612-03



Item No.: 1612-05



Item No.: 1612-06



Item No.: 1612-10

Code	Description	Unit
1612-01	10ml, Ø 33 x 95H mm, ±0.5ml	Each
1612-B	20ml, Ø 42 x 115H mm, ±0.6ml	Each
1612-03	50ml, Ø 49 x 135H mm, ±1.0ml	Each
1612-04	100ml, Ø 58 x 165H mm, ±1.5ml	Each
1612-05	250ml, Ø 80 x 205H mm, ±3.0ml	Each
1612-06	500ml, Ø 104 x 250H mm, ±6.0ml	Each
1612-10	1000ml, Ø 112 x 310H mm, ±10.0ml	Each

PLEASE CONTACT YOUR ACCOUNT MANAGER
TODAY TO PLACE YOUR ORDER