

# Spring Balances

- Available in flat style and tubular
- Graduated in both grams and newtons
- Used to measure the force exerted on an object

Unit: Each

## Flat Style

Code	Description
SP01N	100G/1N
SP10N	1000G/10N
SP20N	2000G/20N
SP25N	2500G/25N
SP50N	5000G/50N
SP100N	10KG/100N
SP100NKG	10KG/100N, Heavy Duty



Code: **SP10N**



Code: **SP50N**

## Tubular Style

Code	Description
SP01NC	100G/1N
SP05NC	500G/5N
SP10NC	1000G/10N
SP25NC	2500G/25N
SP50NC	5000G/50N
SP02NC	250G/2.5N, Aluminium



Code: **SP10NC**



Code: **SP50NC**



**PLEASE CONTACT YOUR ACCOUNT MANAGER  
TODAY TO PLACE YOUR ORDER**

