

Erlenmeyer Conical Flasks

Used to measure, mix and store liquids

- Graduated
- Narrow neck
- Borosilicate glass
- Can be autoclaved or heated over a flame

Unit: Each




Code: FLKEC0025



Code: 1121-0050



Code: 1121-0300



Code: 1121-0500

Code	Description	Unit
FLKEC0025	25ml, Ø 22mm	Each
1121-0050	50ml, Ø 22mm	Each
1121-0300	300ml, Ø 34mm	Each
1121-0500	500ml, Ø 34mm	Each



PLEASE CONTACT YOUR ACCOUNT MANAGER
TODAY TO PLACE YOUR ORDER

