

Separating Funnels

Used in liquid - liquid extractions to separate the components of a mixture

- Borosilicate glass

GLOBAL SHAPE

Code	Description	Unit
1542-A	60ml	Each
1542-A-10	60ml	10/Box
1542-D	125ml	Each
1543-DG	125ml, With Glass Stopcock	Each
1542-E	250ml	Each
1542-E-2	250ml	2/Box
1542-F-2	500ml	2/Box
1542-H	2000ml	Each



Code: 1542-E



Code: 1542-F-2



Item No.: 1543-F

CONICAL SHAPE

Code	Description	Unit
1543-F	500ml, With Stopper and Glass Stopcock	Each

CYLINDRICAL SHAPE

Code	Description	Unit
1543-A	60ml, With Glass Stopcock	Each
1543-D	125ml, With Glass Stopcock	Each