


Pennine

Guedel Airways


Oropharyngeal airways used during anaesthesia or emergencies

- Colour coded bite-blocks for easy size identification
- Smooth airway path for easy cleaning
- Rounded traumatic edges
- Latex free
- Individually wrapped
- Sterile

Unit: Each



Code: **AIRG1P**



Code: **AIRG3P**



Code: **AIRG5P**

Code	Size	Colour
AIRG000P	Size 000, 34mm Length	Light Green
AIRG00P	Size 00, 45mm Length	Blue
AIRG0P	Size 0, 55mm Length	Black
AIRG1P	Size 1, 71mm Length	White
AIRG2P	Size 2, 90mm Length	Green
AIRG3P	Size 3, 100mm Length	Orange
AIRG4P	Size 4, 110mm Length	Red
AIRG5P	Size 5, 115mm Length	Yellow



PLEASE CONTACT YOUR ACCOUNT MANAGER
TODAY TO PLACE YOUR ORDER

