

HAEMOCYTOMETER

A specimen slide which is used to determine the concentration of blood cells in a liquid sample

- Neubauer-improved
- Bright lined with double ruling
- 0.1mm Chamber Depth
- Glass

Unit: Each

Item Code: HCCB



SUPERIOR
MARIENFELD

Haemocytometer Cover Slips

- 20 x 26 x 0.4mm
- Borosilicate 3.3 Glass

Unit: 10/Pack

Item Code: HCS20X26



PLEASE CONTACT YOUR ACCOUNT MANAGER
TODAY TO PLACE YOUR ORDER

