

DRECHSEL BOTTLES AND HEADS

Commonly used to dissolve gasses in a liquid or remove gaseous by-product of a reaction

- Borosilicate glass
- Socket size: 24/29



Item No.: WHFRC200S



Item No.: WHFRC160C

Unit: **Each**

HEADS

Item No.	Length
PLAIN	
DH24-130	130mm, Fits DB24-125
DH24-175	175mm, Fits DB24-250
DH24-215	215mm, Fits DB24-500
SINTERED (<i>Porosity 0</i>)	
DHS24-130-0	130mm, Fits DB24-125
DHS24-175-0	175mm, Fits DB24-250
DHS24-215-0	215mm, Fits DB24-500
SINTERED (<i>Porosity 1</i>)	
DHS24-215-3	215mm, Fits DB24-500

BOTTLES

Item No.	Volume Capacity
DB24-125	125ml
DB24-250	250ml
DB24-500	500ml



**PLEASE CONTACT YOUR ACCOUNT MANAGER
TODAY TO PLACE YOUR ORDER**

