

# MICROSCOPE GLASS SLIDES

Ideal for holding cultures for microscopic examination

- Rectangular piece of glass on which specimens can be mounted for microscopic study
- Thin, flat and resistant to corrosion or fogging

Unit: 50/Pack

## Plain Slides

Item Code	Description
7101-A	Thickness: 0.8-1.0mm, 76.2 x 25.4mm, Pathology Grade, Interleaved
7101-1B	Thickness: 1.0-1.2mm, Pathology Grade, Non-Interleaved
7101-BP	Thickness: 1.0-1.2mm, 76.2 x 25.4mm, Pathology Grade, Interleaved
7103	Thickness: 1.2-1.5mm, with Single Concave Cavity, Interleaved
7104	Thickness: 1.2-1.5mm, with Double Concave Cavity, Interleaved




Item Code: 7101-A



Item Code: 7103



Item Code: 7105-1



Item Code: 7107-PPN

## Frosted Slides

- Frosted one end on both sides
- Thickness: 1.0-1.2mm

Item Code	Description
7105-1	76.2 x 25.4mm, Non-Interleaved
7105-PPA	76.2 x 25.4mm, Pathology Grade, Interleaved
7107-PPN	76.2 x 25.4mm, Pathology Grade, Interleaved
7041-50	76 x 50mm, Double Width



PLEASE CONTACT YOUR ACCOUNT MANAGER  
TODAY TO PLACE YOUR ORDER

