

Honeywell

H910 Plus

N95 Flat Fold Style Respirators

Ideal for Virus Protection, Agriculture, Chemical, Construction, Food Services, Manufacturing, Mining and Utilities

- Approved for 10 times the exposure standard of dusts produced by grinding, drilling, certain applications of brazing and welding etc.
- Shaped for better fit when worn in combination with eyewear
- Low breathing resistance encourages longer wear time by improving the wearer's comfort
- Humidity resistant and hypoallergenic filter media for maximum respirator service life
- Durable latex free suspension straps to eliminate pre-stretching
- No staples, ideal for use in industries such as food
- No restriction of vision or communication



Hazards:

- Airborne particulates - dust/mist/fume

Standards:

- TGA Approved 335988 Medical Device Class 1
- Meets NIOSH 42 CFR 84 N95 requirements for a minimum 95% filtration efficiency against solid and liquid aerosols that do not contain oil
- NIOSH Approved No.: TC-84A-8480

Technical Details:

Attachment Information:	Twin elastic latex free welded straps (staple free)
Outer Shell:	Polypropylene nonwoven
Filter Media:	Polypropylene nonwoven fibres electrostatically loaded making them attractive to particles
Nose Bridge:	Constructed of high strength polypropylene/steel composite increasing bending resistance capacity compared with metal nose clips (contains no aluminium)
Shelf Life:	3 years from date of manufacture

Item Code	Description	Unit
N95R	Head Band with Metal Nose Bridge	50/Box



PLEASE CONTACT YOUR ACCOUNT MANAGER
TODAY TO PLACE YOUR ORDER

