

# N95 and P2 Facemask

- Level 3 barrier protection
- Head band for more secure fit
- Complies with USA NIOSH or European P2 standards
- Filters out solid and liquid particles, excluding particles containing oils
- Mouldable nose piece for comfortable fit
- Non-therapeutic use



Item Code.: **FMN95DK**



Item Code.: **N95R**



Item Code.: **FMN95RNT** and **FMN95R**



Item Code.: **FMN95RVNT**



Item Code.: **MMM8822**



Item Code.: **LIVN95RV**

| Item Code | Respirator Rating  | Standard       | Shape              | Packing | Qty    |
|-----------|--------------------|----------------|--------------------|---------|--------|
| FMN95DK   | N95                | NIOSH 84A-4175 | Duckbill           | Bulk    | 30/Box |
| N95R      | N95                | NIOSH 84A-8480 | Flat Fold Vertical | Bulk    | 50/Box |
| FMN95RNT  | N95                | NIOSH 84A-5411 | Cupped Cone        | Single  | 20/Box |
| FMN95R    | N95                | NIOSH 84A-3323 | Cupped Cone        | Bulk    | 20/Box |
| FMN95RVNT | N95 with Air Valve | NIOSH 84A-5460 | Cupped Cone        | Single  | 10/Box |
| MMM8822   | P2 with Air Valve  | European P2    | Cupped Cone        | Bulk    | 10/Box |
| LIVN95RV  | P2 with Air Valve  | European P2    | Cupped Cone        | Bulk    | 10/Box |



PLEASE CONTACT YOUR ACCOUNT MANAGER  
TODAY TO PLACE YOUR ORDER