

# Storage Bottles

Ideally used in storing samples in liquid or powder form

- Made of high density polythylene (HDPE)
- Highly resistant to acids, alcohols, bases and esters
- Equipped by a strong handle for transportation

Unit: **Each**

## Wide Neck

Code	Description
A11001	5L
A11002	10L
A11003	15L
A11005	25L



Code: A11002

## Narrow Neck

Code	Description
A11021	5L
A11022	10L
A11024	20L
A11025	25L
A11026	50L



Code: A11026

## Narrow Neck with Tap

Code	Description
A11041	5L
A11042	10L
A11043	20L
A11045	25L
A11046	50L



Code: A11041

### TO PLACE AN ORDER,

📞 1300 548 289 @ sales@livingstone.com.au

📄 1300 780 008 🖥️ www.livingstone.com.au