

# Gas Jars

A laboratory glassware used for collecting gas from experiments

- Made of borosilicate glass
- Chemical and heat resistant

Unit: Each

## With Cover

Code	Description
<b>1473-A</b>	50D x 150H mm
<b>1473-D</b>	50D x 200H mm
<b>1473-E</b>	50D x 300H mm



Code: 1473-A

## Without Cover

Code	Description
<b>1473-B</b>	50D x 150H mm
<b>1473-DWC</b>	50D x 200H mm




Code: 1473-B

## Spare Cover

Code	Description
<b>1473-C</b>	Cover for Item No. 1473-A
<b>1473-DC</b>	Cover for Item No. 1473-D



Code: 1473-C



Code: 1473-DC

## TO PLACE AN ORDER,

 1300 548 289   
  sales@livingstone.com.au  
 1300 780 008   
 www.livingstone.com.au



Follow us on **INSTAGRAM**  
@livingstoneinternational

