

# Infrared Thermometers

Accurate temperature readings

## Forehead/Non-Contact Thermometer

- Quick and accurate, 1 Second testing
- Measures body, surface and ambient temperature
- Can take human body temperature with a 3 - 5cm distance from forehead
- Audible alarm if temperature is more than 38°C
- °C or °F temperature display
- Memorise the last 32 Temperature readings
- Digital display screen, 3 Colour backlight LCD
- Automatic power-off (<30 Seconds)
- Long life (Up to 100,000 Readings)



Code: TFIDMRE20LR

## Ear Thermometer

- Safe, hygienic and easy to read
- With 20 Probe covers
- Measures human body temperature through the ear
- Accurate reading
- With LCD screen
- With beeper function for easy operation
- Auto-off function, with memory



Code: THRMET100A

## Ear Thermometer Probe Covers

- Small, plastic shield used to cover the tip of the thermometer
- Prevents contact between the thermometer and the patient's ear



Code: THRMET100A-PCVR

### Code

TFIDMRE20LR

THRMET100A

THRMET100A-PCVR

### Unit

Each

Each

20/Pack

### TO PLACE AN ORDER,

 1300 548 289     sales@livingstone.com.au

 1300 780 008     www.livingstone.com.au



Follow us on **INSTAGRAM**  
@livingstoneinternational

