

INTRAVENOUS (IV) CATHETERS

Ergonomic product design for one-handed insertion

- Colour coded for easier identification of catheter size
- Transparent catheter hub provides confirmation of successful catheter placement when with drawing the needle
- Can be connected to syringe by removing filter cap to expose luer taper end
- Use of hydrophobic membrane filter eliminates blood leakage
- Pyrogen free
- Sterile - Sterilised by ETO (Ethylene Oxide Processing)



CIV1451N



CIV1651N



CIV1851N



CIV2025N



CIV2225N

Listed by the TGA in the Australian Register of Therapeutic Goods – ARTG No. AUST L **326032**

Code	Description	Unit
CIV1451N	14G x 2 Inches, Orange	50 Pieces/Box
CIV1651N	16G x 2 Inches, Grey	50 Pieces/Box
CIV1851N	18G x 2 Inches, Green	50 Pieces/Box
CIV2025N	20G x 2 Inches, Pink	50 Pieces/Box
CIV2225N	22G x 2 Inches, Blue	50 Pieces/Box

PLEASE CONTACT YOUR ACCOUNT MANAGER TODAY TO PLACE YOUR ORDER



Follow us on **INSTAGRAM**
@livingstoneinternational

