

# Burettes

A laboratory apparatus used in quantitative chemical analysis particularly in measuring the volume of a liquid or a gas

- Made of borosilicate glass
- Available in different graduations for accurate measurement

Unit: Each

## Class AS, Teflon Key, Straight Bore Stopcock


| Code      | Description              |
|-----------|--------------------------|
| 435256619 | 10ml, 0.05ml Graduations |
| 435276723 | 25ml, 0.1ml Graduations  |
| 435256830 | 100ml, 0.2ml Graduations |

## Class B, Glass Key, Straight Bore Stopcock

| Code      | Description              |
|-----------|--------------------------|
| 435206619 | 10ml, 0.05ml Graduations |
| 435106725 | 50ml, 0.1ml Graduations  |
| 435206830 | 100ml, 0.2ml Graduations |

## Class B, Automatic, Dr. Schilling, with Schellbach's Stripe

| Code      | Description              |
|-----------|--------------------------|
| 436176619 | 10ml, 0.05ml Graduations |
| 436176723 | 25ml, 0.1ml Graduations  |
| 436176725 | 50ml, 0.1ml Graduations  |



## TO PLACE AN ORDER,

 1300 548 289   
  sales@livingstone.com.au  
 1300 780 008   
 www.livingstone.com.au



Follow us on **INSTAGRAM**  
@livingstoneinternational

