

## Mortars & Pestles

A tool used for pounding and grinding solid substances into fine powder or paste

- A mortar is a vessel where substances are crushed or grounded
- A pestle is a heavy, blunt tool used to grind substances up
- Smooth, rounded bottoms and wide mouths
- Soda lime glass



Code: **1383-CN2**

Code: **1383-B**

### Unit: 1/Box

| Code    | Description |
|---------|-------------|
| 1383-B  | Ø 75mm      |
| 1383-DN | Ø 120mm     |
| 1383-FN | Ø 180mm     |

### Unit: 2/Box

| Code     | Description |
|----------|-------------|
| 1383-AN2 | Ø 60mm      |
| 1383-BN2 | Ø 75mm      |
| 1383-CN2 | Ø 90mm      |
| 1383-DN2 | Ø 120mm     |
| 1383-EN2 | Ø 150mm     |

**PLEASE CONTACT YOUR ACCOUNT MANAGER  
TODAY TO PLACE YOUR ORDER**

