

Conical Flasks

Also known as "Erlenmeyer Flasks", these type of flasks are used in holding, measuring, mixing, heating and boiling liquid chemical samples

- Narrow neck
- Graduated
- Borosilicate glass



Code: FLKEC0150BN-2



Code: FLKEC0250BN-2



Code: FLKEC1000BN-2

Unit: 2/Box

Code	Volume Capacity	Neck Diameter
1121-0050N2	50ml	22mm
FLKEC0150BN-2	150ml	30mm
FLKEC0250BN-2	250ml	32mm
FLKEC0500BN-2	500ml	32mm
FLKEC1000BN-2	1,000ml	43mm

Unit: 1/Box

Code	Volume Capacity	Neck Diameter
1814-0250N	250ml	41mm
1814-0500N	500ml	44mm

PLEASE CONTACT YOUR ACCOUNT MANAGER
TODAY TO PLACE YOUR ORDER



Follow us on INSTAGRAM
@livingstoneinternational

