

# Digital Blood Pressure Monitor

For easy and accurate blood pressure readings

- Measures blood pressure/pulse rate/HB/WHO
- User friendly
- Cuff size 22 - 33cm
- With talking feature
- Auto power-off 3 minutes



## Measures Range:

Pressure: 20 - 280mmHg

Pulse: 40 - 180 per Minute

## Measuring Accuracy:

Pressure: +/- 3mmHg

Pulse: +/- 5%

Unit: Each

Code: **LIVHEM7322**

**PLEASE CONTACT YOUR ACCOUNT MANAGER  
TODAY TO PLACE YOUR ORDER**