

Rectal Catheter Tubes

For safe and effective draining of waste,
performing lavage or administering medicines

- Length: 40cm
- Medical Grade PVC
- With large lateral eyes and rounded tip
- Colour coded for easy identification
- Disposable
- Sterile

Unit: **Each**



Code: CR18PN



Code: CR20PN



Code: CR22PN



Code: CR24PN



Code: CR26PN



Code: CR28PNN

Code	Description
CR18PN	18FG, Ø 6.00mm, Red
CR20PN	20FG, Ø 6.66mm, Yellow
CR22PN	22FG, Ø 7.33mm, Purple
CR24PN	24FG, Ø 8.00mm, Blue
CR26PN	26FG, Ø 8.66mm, White
CR28PNN	28FG, Ø 9.33mm, Grey

**PLEASE CONTACT YOUR ACCOUNT MANAGER
TODAY TO PLACE YOUR ORDER**