

# Plastic Joint Clips

Holds standard taper glassware together


- Made of plastic
- Does not scratch glassware
- Withstands strong concentrated alkalis and dilute acids



Code: KECK14



Code: KECK19



Code: KECK29



Code	Description	Unit
KECK14	14/23, Yellow	Each

## KECK

Code	Description	Unit
KECK10	14/23, Yellow	Each
KECK19	10/19, Mixed Yellow and Green	Each
LSCJC-24	24/29, Green	10/Pack
KECK24	24/29, Green or Black	Each
KECK29	29/32, Orange	Each

PLEASE CONTACT YOUR ACCOUNT MANAGER  
TODAY TO PLACE YOUR ORDER

