

# Tripple Range Ammeters

Measures the flow of electric current in a circuit

- Perfect for students to estimate numerous types of electric currents
- Durable and portable



Code: AMTRTRIR



Code: AMTRTS5MA

Code	Description
AMTRTS1A	1/10/100mA DC, Square
AMTRTRIR	5/50/500mA DC, 4mm Socket, Round
AMTRTS5A	50/500/5000mA, Square
AMTRTS5MA	50mA, 500mA and 5A, Square
AMTRTS500UA	50/500/5000uA, Square

**PLEASE CONTACT YOUR ACCOUNT MANAGER  
TODAY TO PLACE YOUR ORDER**