



ELECTRONIC BABY GROWTH SCALE MEASURE

Accurately and safely measures the height and weight of infants

- Adjustable infant height measurement from 46-80cm with 1 mm graduation
- 8 Memory storage for data review
- Extra-large LCD display screen for easy reading with white backlight
- Auto-off function for power saving
- Toggle between lb/kg
- Hold function for weight determination
- Maximum weight capacity: 25kg/55lbs with 5g graduations
- Includes 4x AAA batteries



Code	Description	Unit
SCLBBYDHW	LCD Screen 75 x 30mm	Each

PLEASE CONTACT YOUR ACCOUNT MANAGER
TODAY TO PLACE YOUR ORDER



Registered NDIS Provider

