



# Pulse Oximeter

Suitable for clinical and personal use

- Fully self-contained and compact
- Equipped with a function switch, through which the oximeter can be automatically powered off when not in use for more than 5 seconds
- With high visibility four-way OLED double colour display
- The rechargeable battery allows for 20 hours of battery life
- Comes with software, USB Cable, A/C charger, rechargeable battery, lanyard and storage container
- Simply attach the lightweight unit to the finger to obtain instant SPO<sub>2</sub>, pulse and pleth wave readings



Unit: **Each**  
Code: **CMS50E**

**PLEASE CONTACT YOUR ACCOUNT MANAGER  
TODAY TO PLACE YOUR ORDER**

