

European Glass Low Form Beakers

Commonly known as Griffin Beaker

- Most commonly used in compound preparation
- With spout - for easy pouring and transferring of samples
- Graduated - helps measure the volume of samples
- Made of borosilicate glass - resistant to thermal shock

Unit: Each

Code	Volume	Diameter	Height
411010005H	5ml	22mm	30mm
411010010H	10ml	26mm	35mm
411010025H	25ml	34mm	50mm
411010050H	50ml	42mm	60mm
417010100H	100ml	50mm	70mm
417010150H	150ml	60mm	80mm
417010250H	250ml	70mm	95mm
417010400H	400ml	80mm	110mm
417010600H	600ml	90mm	125mm
417010800H	800ml	100mm	135mm
417010940H	1L	105mm	135mm
417010950H	2L	130mm	185mm



Code: 411010005H



Code: 417010100H

Unit: 10/Box

Code	Volume	Diameter	Height
411010050HB	50ml	42mm	60mm
417010250HB	250ml	70mm	95mm
417010400HB	400ml	80mm	110mm



Code: 417010250H



Code: 417010400HB

**PLEASE CONTACT YOUR ACCOUNT MANAGER
TODAY TO PLACE YOUR ORDER**