

# Graduated Measuring Cylinders


Also known as mixing cylinder

- Has a narrow cylindrical body shape
- Made of Polypropylene - excellent chemical resistance
- Autoclavable up to 121°C
- Used for accurate and precise measurement of liquid's volume

**Unit: Each**

## Tall Form, Class B

Code	Volume	Outer Ø	Height
A12121	10ml	13.5mm	140mm
A12122	25ml	18mm	195mm
A12123	50ml	25.5mm	199mm
A12124	100ml	30.5mm	249mm
A12125	250ml	41.5mm	315mm
A12126	500ml	55mm	361mm
A12127	1L	66mm	439mm
A12128	2L	84mm	531mm



## Low Form

Code	Volume	Outer Ø	Height
A12181	10ml	16mm	88mm
A12182	25ml	22mm	107mm
A12183	50ml	29mm	143mm
A12184	100ml	34mm	178mm
A12185	250ml	45.5mm	264mm
A12186	500ml	55.5mm	305mm
A12187	1L	70mm	332mm
A12188	2L	92mm	370mm



**PLEASE CONTACT YOUR ACCOUNT MANAGER  
TODAY TO PLACE YOUR ORDER**