


Serological Graduated Pipettes

Commonly known as "Blown-Out Pipettes"


- Designed to accurately measure and transfer a volume of liquid from one container to another
- Especially calibrated to blow-out the last drop of liquid from the tip to deliver the full volume of the pipette
- Sterile, individually wrapped
- Polystyrene
- Accuracy of $\pm 1\%$ for nominal volume




Code: A03101



Code: A03110



Code: A03102



Code: A03105

Item No.	Volume	Division	Colour	Unit
A03101	1ml	1/100mm	Yellow	Each
A03101-1000	1ml	1/100mm	Yellow	1,000/Ctn
A03102	2ml	1/100mm	Green	Each
A03102K	2ml	1/100mm	Green	500/Ctn
A03105	5ml	1/10mm	Blue	Each
A03105K	5ml	1/10mm	Blue	200/Ctn
A03110	10ml	1/10mm	Orange	Each
A03010	10ml	1/10mm	Orange	25/Bag
A03010K	10ml	1/10mm	Orange	500/Ctn
A03110K	10ml	1/10mm	Orange	500/Ctn

**PLEASE CONTACT YOUR ACCOUNT MANAGER
TODAY TO PLACE YOUR ORDER**