

LIVINGSTONE SURGICAL FACE MASK

Level 3 Barrier Protection

- Non-woven 4-ply face mask that protects against dust, dirt and airborne particles
- Triple pleat style for easy breathing
- Flexible nose piece that contours to the face with a foam strip that prevents fogging of eyewear during use
- Tie back style enables tighter fit to suit users
- Bacterial Filtration Efficiency (BFE): $\geq 98\%$
- Level 3 barrier protection
- Meets EN14683:2019 and AC:2019, type IIR standards
- Disposable for single use only
- Latex free



Listed by the TGA in the Australian Register of
Therapeutic Goods – ARTG No. AUST L **388452**

Code	Description	Unit
FM4PAF	Indigo, TGA Listed	50 Pieces/Box

PLEASE CONTACT YOUR ACCOUNT MANAGER
TODAY TO PLACE YOUR ORDER



Follow us on **INSTAGRAM**
[@livingstoneinternational](https://www.instagram.com/livingstoneinternational)

