


# Autoclave Sterilisation Papers

Hygienic storage of bagged instrument and other clinical items, it allows them to be stored in sterilised state ready to be used


- Made of clear film on one side, and medical craft paper on the other side
- Suitable for sterilising devices and equipment
- Has indicator strip to confirm the sterility of objects by changing colour during the process



Code: **PSF050**



Code: **PSF075**



Code: **PSF100**



Code: **PSF150**



Code: **PSF200**

Unit: **Each**

Code	Description
PSF050	50mm x 200m
PSF075	75mm x 200m
PSF075G	75mm x 100m, Gusted
PSF100	100mm x 200m
PSF100G	100mm x 100m, Gusted
PSF150	150mm x 200m
PSF200	200mm x 200m
PSF150G	150mm x 100m, Gusted

Unit: **2 Rolls/Pack**

Code	Description
PSF050-2	50mm x 20m

**PLEASE CONTACT YOUR ACCOUNT MANAGER  
TODAY TO PLACE YOUR ORDER**

