

# Non-Contact Infrared Thermometer

Provides Fast and Hygienic Non-Contact Temperature Measurement

## Specifications:




- Measures Temperature Range: 32 - 42.9°C
- Operates in:
  - Ambient Temperature: 10 - 40°C
  - Humidity:  $\leq 85\%RH$
- Accuracy:  $\pm 0.3^{\circ}C$
- Automatically stores last 32 temperatures for easy recall
- Adjustable mode operation
- Measuring Site: Forehead
- Reference Body Site: Axillary

## Easy Operation:

1. Press On/Scan button, place 3 - 5cm away from forehead, aims lens toward forehead. Press the thermometer's On/Scan button, and the temperature will be displayed instantly
2. Before taking a temperature, make sure any hair, dirt or dust is removed from forehead and forehead is dry. Wait 10 minutes after such removal prior to measuring temperature

**Provides quick results in just 1 second!**  
**5 Years expected service life!**

## Colour Coded Back Light

	38°C - 42.9°C <b>High Fever</b>
	37.4°C - 37.9°C <b>Low Fever</b>
	32°C - 37.3°C <b>Normal Temp.</b>

Code	Unit
THERMOMTALK	Each



**FREE**  
2x 1.5V AA  
Alkaline Batteries



**PLEASE CONTACT YOUR ACCOUNT MANAGER  
TODAY TO PLACE YOUR ORDER**

