

Non-Adherent Dressing

Ideal for open wounds that have light to moderate exudate

- Made of cotton and acrylic bonded to perforated Polyester film
- Can be applied directly to wound
- Non-linting
- Sterile

poren

Thick dressing suitable for skin tears and superficial skin loss wounds

CODE	SIZE	UNIT
LP050X050L5	5 x 5cm	5/Pack
LP050X050	5 x 5cm	100/Box
LP100X100L5	10 x 10cm	5/Pack
LP100X100N	10 x 10cm	50/Box
LP100X200L5	10 x 20cm	5/Pack

**poren
lite**

Thinner dressing suitable for wounds or as a secondary dressing

CODE	SIZE	UNIT
LPL075X050-5	7.5 x 5cm	5/Box
LPL075X050	7.5 x 5cm	100/Box
LPL075X075-5	7.5 x 7.5cm	5/Box
LPL075X075A	7.5 x 7.5cm	80/Box
LPL075X100L5	7.5 x 10cm	5/Box
LPL075X100	7.5 x 10cm	100/Box
LPL075X100-6	7.5 x 10cm	6 Boxes/Ctn, 100/Box
LPL075X200L5	7.5 x 20cm	5/Box

POREN

Code: LP050X050



POREN LITE

Code: LPL075X050



**PLEASE CONTACT YOUR ACCOUNT MANAGER
TODAY TO PLACE YOUR ORDER**

