

# CERTIFIED LABORATORY TEST SIEVES

Designed for particle size testing of chemicals,  
powders, oils, and other fluids

- Ø 200 x 50(H) mm
- Made of stainless steel

Unit: Each

CODE	DESCRIPTION
SIEVE045	45µm
SIEVE063	63µm
SIEVE075	75µm
SIEVE106	106µm
SIEVE150	150µm
SIEVE212	212µm
SIEVE300	300µm
SIEVE425	425µm
SIEVE600	600µm
SIEVE850	850µm
SIEVE2MM	2mm Aperture Wire Mesh
SIEVE4MM	4mm Aperture Wire Mesh



CODE: **SIEVE045**



CODE: **SIEVE2MM**



CODE: **SIEVERCVR**

## Sieve Receiver

CODE	DESCRIPTION
SIEVERCVR	For 200mm Test Sieve

**PLEASE CONTACT YOUR ACCOUNT MANAGER  
TODAY TO PLACE YOUR ORDER**