

DIAGNOSTIC PEN TORCHES

Used as a secondary light source in emergency first aid situations or during first aid treatment, as well as for inspecting ears and throat regions for inflammation and other symptoms.

- 13.8 x 1.3cm
- 2.2 Volts
- 1m Light range
- LED Light bulb
- 8 Lumens light brightness
- Aluminium Alloy frame



Item No.: **TORCHPL3B**



Item No.: **TORCHPL3B**

Item No.	Unit
TORCHPL3B	Black
TORCHPL3DB	Dark Blue

**PLEASE CONTACT YOUR ACCOUNT MANAGER
TODAY TO PLACE YOUR ORDER**