

BOTTLE TOP DISPENSERS

Dispenses chemicals and other liquids precisely from one container to another

- Can withstand 15°C - 40°C
- Vapour pressure up to 500 mbar
- Up to 2.2g/cm³ density
- Polypropylene (PP)

Item No.	Capacity
----------	----------

WLBTP05-05LV	0.5 - 5mL
--------------	-----------

WLBTP1-10LV	1 - 10mL
-------------	----------

WLBTP05-25LV	2.5 - 25mL
--------------	------------



Item No.: **WLBTP05-25LV**

**PLEASE CONTACT YOUR ACCOUNT MANAGER
TODAY TO PLACE YOUR ORDER**