

MICRO TEST TUBES

Used for sample storage, transportation, sample separation, centrifugation and a variety of laboratory applications.

- Ø 10.50 x 40mm, 1.5mL
- Conical with flat cap
- Polypropylene (PP)
- Autoclavable



Item No.

A1003C-500

Unit

500 Pieces/Bag

**PLEASE CONTACT YOUR ACCOUNT MANAGER
TODAY TO PLACE YOUR ORDER**



Follow us on **INSTAGRAM**
[@livingstoneinternational](https://www.instagram.com/livingstoneinternational)

