


# Blood Collection Tubes

Designed to stabilise and preserve blood samples for analysis

- Used for efficient blood extraction
- Made of plastic tube


Code	Additive	Non-Spill Cap	Usage
BCTSST5 & BCTSST	Clot Activator and Serum Separating Gel	Gold	Chemistry, Serology, and Immunology
BCTEDTA6	K2 EDTA	Royal Blue	Heavy Metal, Trace Metals
BD367525 & BD367839 & BCTEDTA2	K3 EDTA	Purple Lavender	Hematology and Blood Bank
BCTEDTA & BCTEDA9	K3 EDTA	Purple Lavender	Hematology and Blood Bank
BCTSER	Silica	Red	Chemistry, Serology, and Immunology
BCTCIT	Sodium Citrate, Na Citrate 3.2%	Light Blue	Coagulation
BCTFLOX	Sodium Fluoride (Glycolysis Inhibitor) & Potassium Oxalate (Anticoagulant)	Grey	Glucose, Blood Alcohol, and Lactate



Code: **BD367525**

Unit: **100 Pieces/Pack**


Code	Description
BCTCIT	2.7mL, 13 x 75mm
BCTEDTA	4mL, 13 x 75mm
BD367839	4mL, 13 x 75mm
BCTFLOX	5mL, 13 x 100mm
BCTSST5	5mL, 13 x 100mm
BCTEDTA6	6mL, 13 x 100mm
BCTSST	8.5mL, 16 x 100mm
BCTEDTA2	10mL, 16 x 100mm
BCTSER	10mL, 16 x 100mm
BD367525	10mL, 16 x 100mm



Code: **BCTSST5**

Unit: **50 Pieces/Pack**

Code	Description
BCTEDA9	9mL, 16 x 100mm



Code: **BCTCIT**

**PLEASE CONTACT YOUR ACCOUNT MANAGER TODAY TO PLACE YOUR ORDER**