

POLYVINYL CHLORIDE (PVC) TUBINGS

Used in different laboratory and vapour transfer applications

- Flexible like rubber
- Can maintain strength and pliancy at low temperature as low as -45°C
- Clear

Unit: Each



Item No.: TUBPL3X1.0MM



Item No.: TUBPL6X1.5MM

1mm Wall Thickness

Item No.	Description	Length
TUBPL3X1.0MM	Inner Ø 3mm, Outer Ø 5mm	1m
TUBPVC05-2	Inner Ø 5mm, Outer Ø 7mm	2m
TUBPVC06-2	Inner Ø 6mm, Outer Ø 8mm	2m

1.5mm Wall Thickness

Item No.	Description	Length
TUBPL3X1.5MM	Inner Ø 3mm, Outer Ø 6mm	1m
TUBPL3X1.5MM-2	Inner Ø 3mm, Outer Ø 6mm	2m
TUBPL3X1.5MMRL	Inner Ø 3mm, Outer Ø 6mm	10m
TUBPL3X1.5MM-10	Inner Ø 3mm, Outer Ø 6mm	10m
TUBPL5X1.5MM	Inner Ø 5mm, Outer Ø 8mm	1m
TUBPL5X1.5MM-2	Inner Ø 5mm, Outer Ø 8mm	2m
TUBPL5X1.5MM-10	Inner Ø 5mm, Outer Ø 8mm	10m
TUBPL6X1.5MM	Inner Ø 6mm, Outer Ø 9mm	1m
TUBPL08X1.5	Inner Ø 8mm, Outer Ø 11mm	1m
TUBPL08X1.5-2	Inner Ø 8mm, Outer Ø 11mm	2m
TUBPL08X1.5-10	Inner Ø 8mm, Outer Ø 11mm	10m
TUBPL10X1.5	Inner Ø 10mm, Outer Ø 13mm	1m
TUBPVC10-2	Inner Ø 10mm, Outer Ø 13mm	2m
TUBPL16X1.5	Inner Ø 16mm, Outer Ø 19mm	1m
TUBPL16X1.5-2	Inner Ø 16mm, Outer Ø 19mm	2m
TUBPL16X1.5-10	Inner Ø 16mm, Outer Ø 19mm	10m
TUBPL16X3.0	Inner Ø 16mm, Outer Ø 22mm	1m
TUBPL19X1.5RL	Inner Ø 19mm, Outer Ø 22mm	10m
TUBPL19X1.5	Inner Ø 20mm, Outer Ø 22mm	1m
TUBPL19X1.5-2	Inner Ø 20mm, Outer Ø 22mm	2m
TUBPL19X1.5-10	Inner Ø 20mm, Outer Ø 22mm	10m
TUBPL25X1.5	Inner Ø 25mm, Outer Ø 28mm	1m
TUBPL25X1.5-2	Inner Ø 25mm, Outer Ø 28mm	2m
TUBPL25X1.5-10	Inner Ø 25mm, Outer Ø 28mm	10m

2mm Wall Thickness

Item No.	Description	Length
TUBPL10X1.5-2	Inner Ø 10mm, Outer Ø 14mm	2m
TUBPL10X1.5-10	Inner Ø 10mm, Outer Ø 14mm	10m
TUBPL12X1.5	Inner Ø 12.5mm, Outer Ø 16.5mm	2m
TUBPL12X1.5-10	Inner Ø 12.5mm, Outer Ø 16.5mm	10m

2.5mm Wall Thickness

Item No.	Description	Length
TUBPL12X2.5	Inner Ø 12mm, Outer Ø 17mm	1m
TUBPL12X2.5-2	Inner Ø 12mm, Outer Ø 17mm	2m
TUBPL12X2.5-10	Inner Ø 12mm, Outer Ø 17mm	10m
TUBPL16X2.5	Inner Ø 16mm, Outer Ø 21mm	1m
TUBPL16X2.5-2	Inner Ø 16mm, Outer Ø 21mm	2m
TUBPL16X2.5-10	Inner Ø 16mm, Outer Ø 21mm	10m
TUBPL19X2.5	Inner Ø 19mm, Outer Ø 24mm	1m
TUBPL19X2.5-2	Inner Ø 19mm, Outer Ø 24mm	2m
TUBPL19X2.5-10	Inner Ø 19mm, Outer Ø 24mm	10m

3mm Wall Thickness

Item No.	Description	Length
TUBPVC08	Inner Ø 8mm, Outer Ø 14mm	1m
TUBPL8X3MM-2	Inner Ø 8mm, Outer Ø 14mm	2m
TUBPVC08-5	Inner Ø 8mm, Outer Ø 14mm	5m
TUBPL8X3MM-10	Inner Ø 8mm, Outer Ø 14mm	10m
TUBPVC08-2	Inner Ø 8mm, Outer Ø 14mm	2m
TUBPVC08-10	Inner Ø 8mm, Outer Ø 14mm	10m
TUBPL16X3.0-2	Inner Ø 16mm, Outer Ø 22mm	2m
TUBPL16X3.0-10	Inner Ø 16mm, Outer Ø 22mm	10m
TUBPL25X3.0	Inner Ø 25mm, Outer Ø 31mm	1m
TUBPL25X3.0-2	Inner Ø 25mm, Outer Ø 31mm	2m
TUBPL25X3.0-10	Inner Ø 25mm, Outer Ø 31mm	10m
TUBPL32X3.0	Inner Ø 32mm, Outer Ø 38mm	1m
TUBPL32X3.0-2	Inner Ø 32mm, Outer Ø 38mm	2m
TUBPL32X3.0-10	Inner Ø 32mm, Outer Ø 38mm	10m

5.5mm Wall Thickness

Item No.	Description	Length
TUBPL40X5.5	Inner Ø 40mm, Outer Ø 51mm	1m
TUBPL40X5.5-2	Inner Ø 40mm, Outer Ø 51mm	2m
TUBPL40X5.5-10	Inner Ø 40mm, Outer Ø 51mm	10m

6mm Wall Thickness

Item No.	Description	Length
TUBPL09X06-2	Inner Ø 9mm, Outer Ø 21mm	2m
TUBPL09X06-10	Inner Ø 9mm, Outer Ø 21mm	10m

**PLEASE CONTACT YOUR ACCOUNT MANAGER
TODAY TO PLACE YOUR ORDER**

