

MULTIPLE REAGENT TEST STRIPS

Used for reliable frontline test and detection of a broad range of conditions, such as urinary tract infections, diabetes and kidney disorders

- The strips can be read visually, requiring no additional laboratory equipment for testing
- A user friendly chart on the back of the bottle makes reading and interpreting results easy

Strips Measure:

Glucose
Bilirubin
Ketone
Nitrite
Blood
pH
Urobilinogen
Protein
Specific gravity
Leucocytes



Item No.	Description	Unit
MULTISTX10LIV	10SG for Urinalysis	100 Strips/Vial

BUY FROM OUR WEBSTORE WWW.LIVINGSTONE.COM.AU
AND GET THE LATEST RANGE AND BEST PRICES