

FILTER FLASKS

Used in vacuum filtration processes to isolate solids from liquid solutions during chemical reactions

- Borosilicate glass
- Remove impurities and particulate matter
- Essential for gravimetric analysis

LIVINGSTONE

Item No.	Capacity	Unit
1431-A	100mL	Each

LINCON

Item No.	Capacity	Unit
1431-BN2	250mL	2 Pieces/Box
074.402.05	500mL	1 Piece/Box

SIMAX®

Item No.	Capacity	Unit
412021100	100mL	Each
412001952	3,000mL	Each
412001950	2,000mL	Each

SUPERIOR MARIENFELD

Item No.	Capacity	Unit
F8.180-24	1,000mL	Each



Item No.: 1431-A



Item No.: 412001952

BUY FROM OUR WEBSTORE WWW.LIVINGSTONE.COM.AU
AND GET THE LATEST RANGE AND BEST PRICES



Follow us on INSTAGRAM
@livingstoneinternational

