

# LINCON

## NITRILE EXAMINATION GLOVES

Provide Superior Resistance to Cuts, Piercing and Blood-Borne Pathogens

- Made from non latex synthetic chemical resistant material without natural rubber latex that minimises allergic reactions to the user
- 4 times stronger than latex or vinyl gloves
- Excellent tactile sensitivity and high level of dexterity for prolonged use
- Features textured fingertips and a beaded cuff for easy donning
- Ambidextrous design
- Conforms to ASTM D6319 and/or HACCP standards
- Cuff Length: 245mm
- Fingertip thickness: Single wall thickness 0.05 mm



### Individual Box

Code	Description	Unit
GLVNT2XS	Extra Small	200/Box
GLVNT2S	Small	200/Box
GLVNT2M	Medium	200/Box
GLVNT2L	Large	200/Box
GLVNT2XL	Extra Large	200/Box

### Carton

Code	Description	Unit
GLVNT2XS-10	Extra Small, 200/Box	10 Boxes/Carton
GLVNT2S-10	Small, 200/Box	10 Boxes/Carton
GLVNT2M-10	Medium, 200/Box	10 Boxes/Carton
GLVNT2L-10	Large, 200/Box	10 Boxes/Carton
GLVNT2XL-10	Extra Large, 200/Box	10 Boxes/Carton

BUY FROM OUR WEBSTORE  
WWW.LIVINGSTONE.COM.AU

AND GET THE LATEST RANGE AND BEST PRICES

