

# LIVINGSTONE AUTOCLAVE/PRESSURE STEAM STERILISERS

Eliminate Microorganisms such as Bacteria and Spores

- Aluminium deep drawn seamless construction
- Sturdy, safe, rust proof, portable, polished shiny finish
- Max pressure 15 PSI/1.5 bar, max temp 121 °C
- Equipped with steam release/control valve, High quality heat resistant bakelite handles
- Pressure gauge with colour coding to display sterilisation zones, 2 2 graduated in kg/cm and lb/inch
- Equipped with safety pressure release weight valve



**ATCLV2100**



**ATCLV910**

Item No.	Description	Unit
ATCLV2100	11L, Ø 9 x 11" Height, Non-Electric, Pressure Cooker Type	Each
ATCLV910	17L, Ø 12 x 10" Height, Electric, 220v/50Hz, Single Phase Wingnut/Knob Type	Each

PLEASE CONTACT YOUR ACCOUNT MANAGER  
TODAY TO PLACE YOUR ORDER



Follow us on **INSTAGRAM**  
[@livingstoneinternational](https://www.instagram.com/livingstoneinternational)

