

LIVINGSTONE CARDIOLOGY LEGEND 3 DUAL HEAD ACOUSTIC STETHOSCOPE

Listen Beyond the Surface

- Acoustic performance rated 8/10, offering noticeably clearer low-frequency sound for detecting crucial cardiac assessments.
- Experience unparalleled convenience with the paediatric side effortlessly transitioning into a traditional open bell by removing the single-piece diaphragm membrane.
- Constructed from robust, acoustically sensitive material, the diaphragm membrane is built to withstand daily use.
- Polished stainless steel finish.
- Included: soft and hard sealing ear tips, a non-chill bell sleeve, and instruction manual.

Removable single-piece diaphragm membrane for use as a traditional open bell



DIA01040 Pediatric and Adult Stainless Steel Stethoscope

Hard-working, dependable, clinical tools for physical assessment and monitoring of your patients, one long shift after another.

| Model | STELLEG3BLK | STELLEG3BLU | STELLEG3BRG | STELLEG3GRY | STELLEG3PNK | STELLEG3RED |
|-----------------------------|-------------------------|-------------------------|-------------------------|-------------------------|-------------------------|-------------------------|
| Colour | Black | Blue | Burgundy | Grey | Pink | Red |
| Approx. Length | 31"/78CM | 31"/78CM | 31"/78CM | 31"/78CM | 31"/78CM | 31"/78CM |
| Approx. Weight | 181g | 181g | 181g | 181g | 181g | 181g |
| Warranty | 2 Years | 2 Years | 2 Years | 2 Years | 2 Years | 2 Years |
| Diaphragm | Turnable (2) | Turnable (2) | Turnable (2) | Turnable (2) | Turnable (2) | Turnable (2) |
| Approx. Chestpiece Diameter | Double-Sided Adult Side | Double-Sided Adult Side | Double-Sided Adult Side | Double-Sided Adult Side | Double-Sided Adult Side | Double-Sided Adult Side |
| | 1.7"/4.4cm | 1.7"/4.4cm | 1.7"/4.4cm | 1.7"/4.4cm | 1.7"/4.4cm | 1.7"/4.4cm |
| | Pediatric Side | Pediatric Side | Pediatric Side | Pediatric Side | Pediatric Side | Pediatric Side |
| | 1.4"/3.5cm | 1.4"/3.5cm | 1.4"/3.5cm | 1.4"/3.5cm | 1.4"/3.5cm | 1.4"/3.5cm |
| Features | Convertible Open Bell | Convertible Open Bell | Convertible Open Bell | Convertible Open Bell | Convertible Open Bell | Convertible Open Bell |

Photos



BUY FROM OUR WEBSTORE
WWW.LIVINGSTONE.COM.AU
AND GET THE LATEST RANGE AND BEST PRICES



Follow us on INSTAGRAM
@livingstoneinternational

