



MEASURING CYLINDERS

Essential laboratory tool designed for the accurate measurement of liquid volumes

- Made of European borosilicate glass
- Tall form
- With hexagonal base and spout

Unit: **Each**

Class A

Blue Stain Graduation

Item No.	Capacity	Dimensions
432110923	25mL	Ø 21.3 x 160(H) mm
432111125N	50mL	Ø 26 x 195(H) mm
432111130	100mL	Ø 31.3 x 240(H) mm
432111546	2L	Ø 85 x 535(H) mm



Item No.: **432110923**

Class B

Amber Stain Graduation

Item No.	Capacity	Dimensions
432150818	5mL	Ø 12.9 x 115(H) mm
432151125	50mL	Ø 26 x 195(H) mm
432151343H	500mL	Ø 53.2 x 380(H) mm
432151546	2L	Ø 85 x 535(H) mm



Item No.: **432150818**

BUY FROM OUR WEBSTORE WWW.LIVINGSTONE.COM.AU
AND GET THE LATEST RANGE AND BEST PRICES



Follow us on **INSTAGRAM**
[@livingstoneinternational](https://www.instagram.com/livingstoneinternational)

