

# SIMAX<sup>®</sup>

## REAGENT BOTTLES

Guarantee safe chemical storage and handling while protecting against contamination and maintaining reagent purity

- Made of European borosilicate glass
- Graduated
- Leak-proof and chemical resistant
- Blue Polypropylene (PP) screw cap and outlet ring

Unit: **Each**



Item No.: **DB0100S**  
**Clear Glass**



Item No.: **DB2000SA**  
**Amber Glass**

### Clear Glass

Item No.	Dimensions	Capacity	Screw Cap
DB0050	Ø 46 x 91(H) mm	50mL	GL32
DB0100S	Ø 56 x 105(H) mm	100mL	GL45
DB0250S	Ø 70 x 143(H) mm	250mL	GL45
DB0500S	Ø 86 x 181(H) mm	500mL	GL45
DB1000S	Ø 101 x 230(H) mm	1L	GL45
DB2000S	Ø 136 x 265(H) mm	2L	GL45

### Amber Glass

Item No.	Dimensions	Capacity	Screw Cap
DB0500SA	Ø 86 x 18(H) mm	500mL	GL45
DB2000SA	Ø 136 x 265(H) mm	2L	GL45

BUY FROM OUR WEBSTORE [WWW.LIVINGSTONE.COM.AU](http://WWW.LIVINGSTONE.COM.AU)  
AND GET THE LATEST RANGE AND BEST PRICES



Follow us on **INSTAGRAM**  
[@livingstoneinternational](https://www.instagram.com/livingstoneinternational)

