

LIVINGSTONE SPECIMEN CONTAINERS

Protect every sample to prevent contamination, ensure accurate identification and safe transport of medical specimens

- Made of Polypropylene (PP)
- With screw cap
- Flat bottom



Item No.: **TP10065S**



Item No.: **TP5744FN**

Labelled

Item No.	Description	Unit
TP5744FN	60 - 70mL Capacity, Brown Screw Cap	550 Pieces/Ctn
TP10065SL	250mL Capacity, Yellow Screw Cap	147 Pieces/Ctn

Unlabelled

Item No.	Description	Unit
TS10844UUN	120mL Capacity, Natural Screw Cap	264 Pieces/Ctn
TP10065UU	250mL Capacity, Natural Screw Cap	147 Pieces/Ctn
TP10065S	250mL Capacity, Yellow Screw Cap	147 Pieces/Ctn
TP12590SU	500mL Capacity, Yellow Screw Cap	72 Pieces/Ctn

BUY FROM OUR WEBSTORE WWW.LIVINGSTONE.COM.AU
AND GET THE LATEST RANGE AND BEST PRICES



Follow us on **INSTAGRAM**
[@livingstoneinternational](https://www.instagram.com/livingstoneinternational)

