

# LIVINGSTONE HYDROMETER JARS

Holds liquid samples for accurate and consistent measurements of chemical substances and beverages, ensuring reliable results

- 1L Capacity
- Base: Ø 101mm
- Inner: Ø 60 x 380(H) mm
- Outer: Ø 65mm
- Made of glass



Item No.	Unit
BACHJ.65	Each

BUY FROM OUR WEBSTORE [WWW.LIVINGSTONE.COM.AU](http://WWW.LIVINGSTONE.COM.AU)  
AND GET THE LATEST RANGE AND BEST PRICES



Follow us on **INSTAGRAM**  
[@livingstoneinternational](https://www.instagram.com/livingstoneinternational)

