



DERMAL FILLER CANNULA

Dome-Shaped Tip Allows Easy Entry For Reducing Bruising

- **Enhanced Precision & Safety:** Patented rounded tip reduces insertion force and lowers the risk of bruising during dermal filler procedures
- **Improved Glide & Comfort:** Ultra-thin wall (UTW) design and proprietary surface coating promote smooth movement and better patient comfort
- **Accurate Delivery:** Clear side-port marking allows for precise control and placement of filler material
- **Minimised Product Waste:** Engineered for efficient filler use, helping reduce wastage and improve procedure outcomes
- **Complete Packaging Solution:** Each cannula is individually packed in a protective capsule and includes a sharp conducting needle for convenience
- **Optimised for Aesthetic Procedures:** Specifically designed to support high-precision dermal filler applications, making it ideal for cosmetic and minimally invasive treatments
- **Individually Packed:** Each sterile pack includes a set of 1 Cannula and 1 Stylet that is ready for immediate use



Item No.	Description	Unit
PRC-22050N	22G x 50mm, Black, With Stylet and Cannula, 1 Set/Pack	20 Packs/Box
PRC-22070N	22G x 70mm, Black, With Stylet and Cannula, 1 Set/Pack	20 Packs/Box
PRC-25038N	25G x 38mm, Orange, With Stylet and Cannula, 1 Set/Pack	20 Packs/Box
PRC-25050N	25G x 50mm, Orange, With Stylet and Cannula, 1 Set/Pack	20 Packs/Box
PRC-27038N	27G x 38mm, Grey, With Stylet and Cannula, 1 Set/Pack	20 Packs/Box
PRC-27050N	27G x 50mm, Grey, With Stylet and Cannula, 1 Set/Pack	20 Packs/Box

BUY FROM OUR WEBSTORE
WWW.LIVINGSTONE.COM.AU
AND GET THE LATEST RANGE AND BEST PRICES

