


LIVINGSTONE **AUTOCLAVE TAPES**

Designed for sealing autoclave bags containing instruments and samples that will undergo sterilisation

With Indicator

- Changes colour to indicate successful sterilisation
- Activates after 20 minutes of steam sterilisation at 121 °C – 125 °C
- Line turns black after steam sterilisation
- "GAS" text turns green after gas sterilisation
- Class 4 chemical indicator

Item No.	Description	Unit
TA12X30	12mm x 30m	Each
TA12X30-10	12mm x 30m	10 Pcs/Box
TA12X50	12mm x 50m	Each
TA12X50-10	12mm x 50m	10 Pcs/Box
TA19X30	19mm x 30m	Each
TA19X30-10	19mm x 30m	10 Pcs/Box
TA19X50	19mm x 50m	Each
TA19X50-10	19mm x 50m	10 Pcs/Box
TA19X55	19mm x 55m	Each
TA25X30	25mm x 30m	Each
TA25X30-10	25mm x 30m	10 Pcs/Box
TA25X50	25mm x 50m	Each
TA25X55	25mm x 55m	Each
TA25X55-10	25mm x 55m	10 Pcs/Box




Item No.: **TA12X30**

Without Indicator

- Good for general lab use, industrial parts, or items not requiring confirmed sterility
- Can be used when internal chemical or biological indicators are already in place for efficient and cost-conscious process

Item No.	Description	Unit
TA12X30NI	12mm x 30m	Each
TA25X30NI	25mm x 30m	Each
TA25X55NI	25mm x 55m	Each



Item No.: **TA25X55NI**

**BUY FROM OUR WEBSTORE WWW.LIVINGSTONE.COM.AU
 AND GET THE LATEST RANGE AND BEST PRICES**

