

LIVINGSTONE MULTI-SAMPLE BLOOD COLLECTION NEEDLES


Enables multiple blood samples from a single puncture and reduces patient discomfort

- Provides instant confirmation of needle placement through a unique visual flashback indication when attached to a vacutainer
- Sterile double-ended needle with a sleeve-covered safety valve
- Advanced low-angle bevel design and silicone lubricant coating ensure gentle and smooth vein entry
- Draws blood directly from the vein into vacuum collection tubes to prevent exposure
- Can be used with multiple collection tubes after venipuncture
- Pen type
- Individually packed

Unit: **100 Pieces/Box**



Item No.: **LBCN20X1**



Item No.: **LBCN21X15**



Item No.: **LBCN22X1**

Item No.	Description	Colour
LBCN20X1	20G x 1" (25mm)	Yellow
LBCN20X15	20G x 1.5" (38mm)	Yellow
LBCN21X1	21G x 1" (25mm)	Green
LBCN21X15	21G x 1.5" (38mm)	Green
LBCN22X1	22G x 1" (25mm)	Black
LBCN22X15	22G x 1.5" (38mm)	Black

BUY FROM OUR WEBSTORE WWW.LIVINGSTONE.COM.AU
AND GET THE LATEST RANGE AND BEST PRICES



Follow us on **INSTAGRAM**
@livingstoneinternational

