

LIVINGSTONE **BRIX REFRACTOMETERS**

Convenient tool for quick measurement of concentration of substance in liquid and semi-solid samples

- 175g, Lightweight
- No power needed
- Result in seconds



Item No.: **BRXREF104**



Item No.: **BRXREF108**

Brix

Item No.	Size	Range	Accuracy
BRXREF104	Ø 29 x 148(L) mm	28 – 62%	0.20%
BRXREF105	Ø 29 x 148(L) mm	40 – 82%	0.50%
BRXREF106	Ø 29 x 148(L) mm	58 – 92%	0.50%
BRXREF109	Ø 29 x 148(L) mm	0 – 90%	0.50%
BRXREF101	Ø 29 x 160(L) mm	0 – 10%	0.10%
BRXREF102	Ø 29 x 160(L) mm	0 – 20%	0.10%
BRXREF103	Ø 29 x 160(L) mm	0 – 32%	0.20%

Brix with ATC (Automatic Temperature Compensation)

Item No.	Size	Range	Accuracy
BRXREF108	Ø 29 x 148(L) mm	0 – 80%	0.50%
BRXREF114	Ø 29 x 148(L) mm	8 – 62%	0.20%
BRXREF111	Ø 29 x 160(L) mm	0 – 10%	0.10%
BRXREF113	Ø 29 x 160(L) mm	0 – 32%	0.20%

**BUY FROM OUR WEBSTORE WWW.LIVINGSTONE.COM.AU
AND GET THE LATEST RANGE AND BEST PRICES**



Follow us on **INSTAGRAM**
[@livingstoneinternational](https://www.instagram.com/livingstoneinternational)

