



LABORATORY CLING FILMS

Used to seal lab containers, preventing sample evaporation, and contamination

- Plasticisers free
- Stretches up to 200%
- Clings to irregular shapes and surfaces
- Polar substance resistance up to 48 hours
- Made primarily from polyolefins and paraffin waxes

Unit: **Each**



Item No.: **PM996**



Item No.: **PM992**

Item No.	Size
PM996	10cm x 38m
PM992	5cm x 76m
PM992S	5cm x 76m (slight imperfection)

BUY FROM OUR WEBSTORE WWW.LIVINGSTONE.COM.AU
AND GET THE LATEST RANGE AND BEST PRICES



Registered NDIS Provider

